

Scope of P2106-RT

Attention:

It is essential to apply a hydrolysis step for this ring test. An analysis of the esters and the free acids only results in an underestimation of the spiked levels! The test material MUST be analysed twice - with and without hydrolysis.

Free acids, esters and conjugates of a selection of the acidic herbicides mentioned below are spiked to the test material of P2106-RT.

4-CPA
2,4-D
2,4-DB
2,4-Dichlorprop
2,4,5-T
2,4,5-TP
Clopyralid
Dicamba
Fenoxaprop-P
Fluazifop
Fluroxypyr
Haloxypop
MCPA
MCPB
Triclopyr